



Reduce cost with Invacare high flow concentrator

The *Invacare Platinum 9* offers providers a value added and reliable global solution to treat patients with high flow requirements. Use as a stand alone device or couple together with *HomeFill* for complete patient care.

While optimising oxygen concentration, *Platinum 9* minimises costs associated to maintenance and delivery. It allows patients to stay at home comfortably and safely thanks to a series of integrated alarms. Users will appreciate its simplicity of use combined with maximum security.



Trusted Solutions, Passionate People



Yes, you can.®



High performance and reliability combined with simplicity

Powerful and reliable, *Invacare Platinum 9* is also user friendly. With its 9 litre oxygen capacity and 96-87% purity, **Platinum 9** is dedicated to patients with high oxygen needs. Easy to manoeuvre thanks to its integrated handle and casters, **Platinum 9** is also easy to use and guarantees extra safety for its user.

Value added associated to cost control

Cost linked to inventory is kept to a minimum thanks to the small number of spare parts and specific accessories. With no maintenance for three years, Invacare offers you the best value added concentrator.

Features and accessories



Secured power cord connection



Only one filter to change when dirty



HomeFill connection port

Technical data

Platinum 9	365 mm	467 mm	24 kg	230V, 50Hz	475 W
Platinum 9	1 - 9 l/min	62 kPa	< 47 dB(A)	87-96%	Up to 1230m ISO8359, MDD

SensO₂-System

O₂-concentration

higher than 85%
- green light -
Normal

between 73% and 85%
- yellow light -
Check filter and flow rate

lower than 73%
- red light & audible alarm -
Call service technician

Service Interval

No provider routine preventative maintenance for 3 years, replace filters as required.

Imported & Distributed by:



SOUTHSIDE BIOMEDICAL SALES AND SERVICES CO. LTD.

**#903 9TH FLOOR PHILFLEX BAY CENTER,
CORAL WAY DRIVE, PASAY CITY**

8-461-3468/ 0917 728 0229

info@southsidebiomedical.com