

AM-2000H 12.1 inch Patient monitor



Feature

- 1.Suitable for Bedside,ICU,Operation room
- 2.Suitable for adult, pediatric and neonatal patients
- 3.Support various modules including 3/5/12 Leads ECG,RESP,SPO2, NIBP,IBP,CO2,AG etc.
- 4.Multi-display:Large font, Standard,Full-ledas ,trend,oxyCRG.
- 5.Built-in detachable rechargeable lithium battery(about 4h)
- 6.Support multi-language,Central Monitoring system
- 7.Unique voice and indicator alarm
- 8.Maximum 720 hour graphic and tabular trends of all parameters
- 9.Arrhythmia analysis, ST calculation, Pace analysis
- 10.Anti-ESU, anti- defibrillator

Options

- .Touchscreen
- .Printer
- .Dual IBP
- .ETCO2
- .Anesthesia multi-gas
- .CO(Cardiac output)
- .Central monitor system
- .WIFI
- .12 Lead ECG
- .Wall mounted basket
- .Monitor trolley

