



DAYPRO TECH



Catalogue for Oxygen Concentrator ZY-1S&1Z&1J

► Catalogue for Oxygen Concentrator ZY-1S&1Z&1J

Established in 2013, Beijing Day-Pro Technology Co., Ltd. is a leading professional company in the field of laboratory and gas generating equipment.

Pursing high quality and excellent performance is our motto.

Our advantages are as the following:

- We have very strong technical team of 150 staff.
- We have our own factory to design and produce according to your requirements.
- We have strong team to manufacture and service.
- We provide customized and complete solution.
- We focus on long-term business cooperation with our client !



ZY-1S

ZY-1J

ZY-1Z


Technical Parameter

Technical Parameter	ZY-1S	ZY-1S (Upgrade)	ZY-1Z	ZY-1Z (Upgrade)	ZY-1J	ZY-1J (Upgrade)
Gross weight	7.1KG	7.1KG	7KG	7KG	7KG	7KG
Net weight	6.5KG	6.5KG	6.5KG	6.5KG	6.5KG	6.5KG
Product Size	266*230*295 (mm)	266*230*295 (mm)	284*187*302 (mm)	284*187*302 (mm)	280*192*300 (mm)	280*192*300 (mm)
Box gauge	370*300*375 (mm)	370*300*375 (mm)	365*270*365 (mm)	365*270*365 (mm)	365*270*400 (mm)	365*270*400 (mm)
Power	220V/110V	220V/110V	220V/110V	220V/110V	220V/110V	220V/110V
Capacity	240VA	240VA	240VA	240VA	240VA	240VA
Oxygen Purity	90%±3%	90%±3%	90%±3%	90%±3%	90%±3%	90%±3%
Oxygen Flow	1-7L	1-7L	1-7L	1-7L	1-7L	1-7L
Atomization	N/A	Available	N/A	Available	N/A	Available

Features are as the following

- Counts: It can account the time you used.
- Intelligent Voice , It can Voice announce the operation.
- Multiple muffler: It has several noise reducing device which will resulted in very low noise.
- Purity monitoring : You can monitor the oxygen purity and make adjustment.
- Remote Control: You can control and operate the machine far by using the remote controller.

► Working Principle (PSA)



► Warranty : One Year

► Certification: CE

Imported & Distributed by:

SOUTHSIDE
Biomedical

SOUTHSIDE BIOMEDICAL SALES AND SERVICES CO. LTD.

📍 #903 9TH FLOOR PHILFLEX BAY CENTER,
CORAL WAY DRIVE, PASAY CITY

☎ # 8-461-3468 / 0917 728 0229

✉ info@southsidebiomedical.com